









ECLAIRAGE OPERATOIRE ALM ANGENIEUX

ALM ANGENIEUX	PRISMATIC ECA002		PRISMALIX LHX001	
				
<i>Equivalence</i>	<i>SYLVANIA</i>	<i>LAMPE MAQUET</i>		<i>LAMPE MAQUET</i>
Référence	727097	727ECA	727095	727LHX
Culot	ailettes		ailettes	
Volts	24		23	
Watts	120		100	

ALM ANGENIEUX	AXCEL / X-TEN AX4 / AX10		AX8	
				
<i>Equivalence</i>	<i>AX 186762</i>	<i>AX 163508</i>	<i>AX 76554</i>	<i>sablée 39230</i>
Référence	715026	715025	715031	727032
Culot	G6,35	G6,35	G6,35	G6,35
Volts	24	24	24	24
Watts	100	150	250	250

ANGENIEUX	AX 4		AX 14	
				
<i>Equivalence</i>	<i>AX 153092</i>	<i>sablée AX 104316</i>	<i>AX153091</i>	<i>sablée AX150394</i>
Référence	727090	727092	727093	727094
Culot	ailettes X511	ailettes X511	ailettes X511	ailettes X511
Volts	24	24	24	24
Watts	70	70	100	100

- VOIR EGALEMENT NOS SUPPORTS DE LAMPES -