






ECLAIRAGE OPERATOIRE Dr MACH

Dr MACH	Praxila / Triaflex / 120F		Soloflex / Uniflex / Trigenflex / Quintaflex 130F/380F/400F/500F/700F	
	<i>Equivalence</i>			
	<i>1240/91</i>	<i>Fil. Axial</i>	<i>1660/91</i>	<i>Fil. Axial</i>
Référence	727002	727003	727004	727005
Culot	G6,35	G6,35	G6,35	G6,35
Volts	22,8	22,8	22,8	22,8
Watts	40	40	50	50

Dr MACH	M3	M2 / M5	Praxila Triaflex	Soloflex
				
<i>Equivalence</i>			<i>1210/91</i>	<i>1150/91</i>
Référence	727001	727016	727125	727127
Culot	G6,35	G6,35	Ba15d	Ba20d
Volts	22,8	22,8	24	24
Watts	40	80	40	60