






DENTAIRE ECLAIRAGE OPERATOIRE

				
Référence	727012	727023	727024	664382
Culot	Pk22s	Pk22s	Pk22s	Pk30d
Volts / Ampères	12	24	24	6,6A
Watts	55	65	70	200
Durée de vie (h)	225	-	225	1000
Réf. OSRAM	64151	-	64156	64382

			
Référence	664340	664380	623710
Culot		R7s	R7s
Ampères / Volts		6,6A	25V
Watts	100	200	150
Durée de vie (h)		1000	3000
Réf. OSRAM	64340	64380	-
Réf. GE	-	-	GE Q150T4

			
Référence	614623	664641	614623P
Culot	G6,35	G6,35	GY9,5
Volts	17	24	17
Watts	95	150	95
Durée de vie (h)	2000	2000	2000
Réf. OSRAM	-	64641	-
Réf. PHILIPS	14623	-	14623P