
























ECLAIRAGE OPERATOIRE HANAULUX

HANAULUX	Hospital			Hambourg			Amsterdam
	 	 	 				
<i>Equivalence</i>	16191	16164	18550	18366			
<i>Référence</i>	727110	727112	727125	727015			
Culot	Ba15d	Ba15d	Ba15d	R7s			
Volts	12	24	24	22,8			
Watts	25	25	40	75			

HANAULUX	Palma	Boston, Londres, Vienne			Universal
	 	 	 	 	
<i>Equivalence</i>	18306	16678	18251	dépolie 18250	
<i>Référence</i>	727080	727082	727085	727086	
Culot	Bx22d	Bx22d	Bx22d	Bx22d	
Volts	24	24	220	220	
Watts	40	50	50	50	

	Melbourne Oslo Heidelberg	Londres Francfort	Hanau 2000	HLX 3000
HANAULUX				
<i>Equivalence</i>	18769	18566	15282	XIR001
Référence	727002	727004	727005	727001
Culot	G6,35	G6,35	G6,35	G6,35
Volts	22,8	22,8	22,8	22,8
Watts	40	50	50	40

	Blue 30	Blue 80	Blue 90/130	Blue 100
HANAULUX				
<i>Equivalence</i>	53010	53026	53198	53296
Référence	727174	727175	727176	727177
Culot	cablé	spécial	spécial	spécial
Volts	24	21,5	22,8	22,8
Watts	50	130	90	113