










## ENDOSCOPIE XBO

	Non câblée	Connecteur 5 voies	Cosses anneaux
<b>XBO</b>			
	<i>R100/45</i>	<i>R100/45C</i>	<i>R110/45C<sup>(1)</sup></i>
<b>Référence</b>	<b>510045NC*</b>	<b>510045*</b>	<b>510145*</b>
<b>Volts</b>	<b>14</b>	<b>14</b>	<b>14</b>
<b>Watts</b>	<b>100</b>	<b>100</b>	<b>100</b>
<b>Durée (H)</b>	<b>500</b>	<b>500</b>	<b>500</b>

	Non câblée	Connecteur 5 voies	Cosses anneaux
<b>XBO</b>			
	<i>R180/45</i>	<i>R180/45C</i>	<i>R181/45C<sup>(1)</sup></i>
<b>Référence</b>	<b>518045NC*</b>	<b>518045*</b>	<b>518145*</b>
<b>Volts</b>	<b>14</b>	<b>14</b>	<b>14</b>
<b>Watts</b>	<b>180</b>	<b>180</b>	<b>180</b>
<b>Durée (H)</b>	<b>500</b>	<b>500</b>	<b>500</b>

	Cosses plates		
<b>XBO</b>			<b>XBO</b>
	<b>droite</b>	<b>coudée</b>	
	<i>R181/45C<sup>(1)</sup></i>		<i>R300/60C</i>
<b>Référence</b>	<b>518145CD*</b>	<b>518145CC*</b>	<b>530060*</b>
<b>Volts</b>	<b>14</b>		<b>17,5</b>
<b>Watts</b>	<b>180</b>		<b>300</b>
<b>Durée (H)</b>	<b>500</b>		<b>500</b>

\*Soumis à Eco-contribution

(1) Préciser la longueur des câbles de connexion.