




DENTAIRE

LAMPE A POLYMERISER

LAMPE A POLYMERISER						
	Référence	613865	664617	Réfect. bleu 664613	664624	EJM 615995
Culot	G5,3	G5,3-4,8	G5,3-4,8	G5,3-4,8	G5,3-4,8	GX5,3
Volts	12	12	12	12	12	21
Watts	75	75	75	100	100	150
Heures	50	25	25	25	25	40
<i>Réf. OSRAM</i>	-	64617	64613	64624	64624	-
<i>Réf. PHILIPS</i>	13865	-	-	-	-	5995

LAMPE A POLYMERISER					
	Référence	613294	613298	613528	613165
Culot	GZ4	GZ4	GZ4	GZ4	GZ4
Volts	6	10	6	6	14
Watts	35	52	15	15	35
Heures	35	25	750	750	50
<i>Réf. PHILIPS</i>	13294	13298	13528	13528	13165