







ENDOSCOPIE LAMPES HALOGENES

Largeur x hauteur	10 x 40			
				
Référence	715017	715021	715019	715023
Code ansi	BRL	FCR	EVA	BRJ/EVB
Culot	G6,35	GY6,35	GY6,35	G6,35
Volts	12	12	12	15
Watts	50	100	100	150
Heures	50	50	2000	50
Réf. OSRAM	64610	64625	64623	64633
Réf. Philips	7027	7023	7724	6550

Largeur x hauteur	13 x 45		13 x 50	
				
Référence	715027	715025	715033	715031
Code ansi	FCS	FDV	EHJ	EVC
Culot	G6,35	G6,35	G6,35	G6,35
Volts	24	24	24	24
Watts	150	150	250	250
Heures	50	300	50	300
Réf. OSRAM	64640	64642	64655	64657
Réf. Philips	7158XHP	6680	7748XHP	6958

				
Référence	715080	715090	715100	715110
Code ansi	EFP	EFR	EKE	ELC
Culot	GZ6,35	GZ6,35	GX5,3	GX5,3
Volts	12	15	21	24
Watts	100	150	150	250
Heures	50	50	200	50
<i>Réf. OSRAM</i>	<i>64627</i>	<i>64634</i>	<i>93638</i>	<i>64653</i>
<i>Réf. Philips</i>	<i>6834</i>	<i>6423</i>	<i>13629</i>	<i>GE 37462</i>

FUJINON		
	LMP-HA75	LMP-SD
Référence	715079	715093
Culot	GZ4	GZ6,35
Volts	12	15
Watts	75	150

					
Référence	715042	715041	715140	715115	715120
Code ansi	FDS/DZE	FDT	DNE	DNF	ELV
Culot	GY9,5	GY9,5	GX7,9	GX7,9	GX7,9
Volts-	24	12	120	21	21,5
Watts	150	100	150	150	150
Heures	100	50	15	25	-
Réf. OSRAM	64643	64628	54409	93631	54765
Autres Réf.	-	PHI 5973	SYL 60742	-	SYL 60936