







## LAMPES POUR MICROSCOPES HBO

<b>HBO</b>						
<b>Référence</b>	<b>925132*</b>	<b>925130L1*</b>	<b>925130L2*</b>	<b>925134*</b>	<b>925136*</b>	<b>925138*</b>
Volts	23	39-45	34-39	20	23	61
Watts	50	50	50	100	100	200
Heures	200	100	100	200	300	200
<i>Réf. OSRAM</i>	<i>HBO 50W/3</i>	<i>HBO 50W/ACLI</i>	<i>HBO 50W/ACL2</i>	<i>HBO 100W/2</i>	<i>HBO 103W/2</i>	<i>HBO 200W/4</i>

<b>HXP</b>	
<b>Référence</b>	<b>925140*</b>
Volts	75
Watts	120
Heures	2000
<i>Réf. OSRAM</i>	<i>HXP R120W / 45C VIS</i>

<b>STARTER SPECIAL HBO</b>	
<b>Référence</b>	<b>300192</b>
<i>Réf. OSRAM</i>	<i>ST192</i>

\*Soumis à Eco-contribution