











MICROSCOPES ZEISS

ZEISS	 <p style="text-align: center;">PG22</p>				
Référence	915046	915045	915047	915048	915049
Code ansi / lif	-	-	M185	-	EHE/A145
Culot	PG22	PG22	PG22	PG22	PG22
Volts	6V	6V	12V	12V	12V
Watts	10W	20W	30W	50W	100W
Heures	300	300	50	50	50
Réf. OSRAM	64222	64251	64260	64609	64626
Equivalence	380120-7030D	38-13-87	380120-7040	-	380075-1020

ZEISS					
Référence	925108		925110		925112
Culot	P47d		P47d		
Volts	6V		6V	6V	
Watts	25W		30W	50W	
Heures	80		80		
Equivalence	39-01-53		39-01-58	39-01-86	

ZEISS					9 x 30  <i>Axioscope</i>
	Référence	925102	925104	918023	715006N
Culot	P22d	P25d4	E10	G4	
Volts	6V	6V	6V	6V	
Watts	15W	15W	15W	20W	
Heures	80	80	200	100	
<i>Equivalence</i>	<i>380018-1730</i>	<i>380018-1740</i>	<i>1521</i>	<i>380079-9690</i>	

ZEISS			
	Référence	925038	925050
Culot	Ba15s	Ba15d	Ba20d
Volts	6V	6V	12V
Watts	10W	15W	60W
Heures	600	1000	-
<i>Equivalence</i>	<i>OSRAM 8014</i>	<i>380029-9020</i>	<i>38-02-16</i>