


















MICROSCOPES / COLPOSCOPES



	9 x 38	9 x 30		9 x 30
				
	Préciser Osram-Philips	Filament vertical		
Référence	715002	715006	715009	715008
Culot	G4	G4	G4	G4
Volts	6V	6V	6V	12V
Watts	10W	20W	35W	20W
Heures	300	100	50	2000
<i>Réf. OSRAM</i>	64223	64250	64275	64258
<i>Réf. Philips</i>	5972	7388	-	-


	11,5 x 44		
			
Référence	715017	715019	715021
Culot	G6,35	GY6,35	GY6,35
Volts	12V	12V	12V
Watts	50W	100W	100W
Heures	50	2000	50
<i>Réf. OSRAM</i>	64610	64623	64625
<i>Réf. Philips</i>	7027	7724	7023

NIKON					
		Filament plat			
Référence	925106	925036	925048	925052	715057
Culot	P25d	Ba15d	Ba15d	Ba15d	GZ4
Volts	6V	6V	6V	6V	6V
Watts	15W	15W	15W	30W	20W
Heures	80	-	100	200	50
<i>Equivalence</i>	77903	77910	76906	77904	20533

CARL ZEISS IENA				
Référence	925115	925122	925125	925126
Culot	Z16	PY16-1,25		
Volts	6V	6V	6V	6V
Watts	15W	20W	25W	30W
Heures	100	100		
<i>Equivalence</i>		-	41780-9001	-

WILD			 
Référence	925114	925113	925109
Culot	support spécial		Ba15d
Volts	12V	12V	6V
Watts	30W	50W	20W
Heures	50	50	-
<i>Equivalence</i>	384643		177160

GAMBS	 
Référence	925107
Culot	P45d
Volts	6V
Watts	15W
Heures	100
<i>Equivalence</i>	Gambs

	
Référence	715082
Culot	GZ6,35
Volts	12
Watts	100
Heures	1500
<i>Equivalence Osram</i>	64637